



SafeBlue insemination catheter



SafeBlue hygiene concept

Insemination success is based on three pillars: sperm quality, the moment of AI and, last but not least, hygiene. This is where the **unique SafeBlue concept for hygienic porcine insemination** comes in.

SafeBlue catheters are single packed, sterile and pre-lubricated. The sheath acts as protection foil during storage and especially during the passage through the vulva. The catheter remains free of contamination until insemination.

A number of scientific evaluations and field reports emphasize the positive effect of the SafeBlue concept: **30-100 more piglets per 100 inseminated breeding sows** have been consistently achieved in practice.

SafeBlue catheters are available in several versions:

SafeBlue ClearGlide® traditional form	
Safelock, 100/bag	17106/3001
with handle, 100/bag	17106/4001
SafeBlue ClearGlide® with sealing flange [1]	
Safelock, 100/bag	17106/3005
with handle, 100/bag	17106/4005
SafeBlue Spirette® with handle, 100/bag [2]	
	17102/9561
SafeBlue Foamtip® Safelock, 100/bag [3]	
	17106/9076
ProFit SafeBlue Foamtip® with handle, 1/bag [3]	
	17118/1125
ProFit SafeBlue SoftGilt [4]	
1/bag	17119/1025
with handle, 1/bag	17119/1125
SafeBlue Foamtip®, with extension, 100/bag [5]	
	17107/0133
Flexible extension, 1/bag [6]	
	17108/1001



Minitube extensions...

... are suitable for all insemination catheters. They are used to hang semen tubes on a wire, providing a hygienic connection between the semen dose and the catheter. The flexible extensions can be used for all insemination techniques as the tube/bag can be fixed in almost any position.






(+ Your benefits

- + Very clean insemination with little effort: the sterilized catheter remains clean and free of contamination during storage until insemination
- + The protective sanitary sheath keeps the catheter clean during insertion into the vulva, after which the catheter tip is pushed through the sheath
- + No more fumbling with lubricants
- + Easy to prepare and quick to use: more inseminations per hour per technician
- + Improved hygiene: increased pregnancy rates and litter sizes (30-100 more piglets per 100 inseminated breeding sows)

(★ Application

- Clean the outside of the sow's vulva with a dry paper towel
- Introduce the catheter with the sanitary sheath 5-10 cm; then push the catheter tip through the sheath
- Pull the sheath backwards
- Unscrew the tip of the QuickTip® tube and lock it into the catheter shaft with a quarter turn
- Start insemination

Improved pregnancy rates with AI using SafeBlue catheters



141 more piglets per 100 sows inseminated
 ⇒ 1.41 more piglets per sow (+15.8 %)

In the field trial by Hoy et al. (2005a), the use of SafeBlue resulted in 141 more piglets per 100 sows compared to a conventional insemination catheter.

Literature

Hoy, S., De Alba, C. (2005): Minimizing the risk of contamination during artificial insemination improves reproductive results in sows. *Reproduction in Domestic Animals* Vol. 40, No. 4, pp 390.

Hoy, S., De Alba, C. (2005): Improvement in breeding results by technician skills and hygiene during routine artificial insemination of sows. *Proceedings of the 7th International Conference on Pig Reproduction*, Kerkrade, The Netherlands, pp 134.

Better hygiene during AI improves fertility (2004). *International Pig Topics*, Volume 19, No. 6.

Heinze, A. (2006): „Mehr Aufmerksamkeit für Besamungs-hygiene.“ *Thüringer Landesanstalt für Landwirtschaft*

Hoy, S. (2005a): Die „sterile“ Besamung – Eingeschweißter Besamungskatheter erhöht die Abferkelrate. *BLW 2*.

Hoy, S. (2005b): More piglets with a single wrapped sterilized catheter. *Pig Process*, Volume 21, No. 8.

TA Temmen, F. (2003): Erfahrungsbericht zum Besamungs-katheter „SafeBlue“. *BES Golzow, Deutschland*.

Reicks DL (2003): Bacterial contamination and semen quality. *Proc. Allen D. Lemman Swine Conf.*, 169-170.

Thompson, R. (2000): Transportation, cleaning and disinfection swine health fact sheet. Vol. 2, n°2, January, NPPC.

De Winter PPJ, Verdonck, M., de Krief, A., Devriese LA, Haesebrouck (1992): Endometritis and vaginal discharge in the sow. *Anim Reprod Sci*; 28: 51-58.