



Technical Data Sheet

Collection dummy for boars

Boar dummy height adjustable with gas pressure spring, designed for hygienic and efficient semen collection with maximum comfort for the boar. Heavy steel body coated with robust and removable PUR plastic or stainless steel to provide a hygienic surface, removable mounting element, footrests and base with locking pin.

The boar dummy is constructed so that the plastic cover and mounting can be removed with little effort. This allows for easy cleaning and replacement. BoarMatic mechanics can be added.



Dimensions (L x W x H)

- Dummy 11 inch: 860 x 285/500 x 532-703 mm
- Dummy 9 inch: 860 x 230/450 x 532-703 mm

Weight

- Dummy 11 inch: 80 kg
- Dummy 9 inch: 76 kg

Boar dummy	
11" for mature boars, plastic cover	11100/0207
11" for mature boars, stainless steel body	11100/0210
9" for mature boars, plastic cover	11100/0107
9" for mature boars, stainless steel body	11100/0108



The mounting element can be easily removed in 2 steps



Accessories and replacement parts

Mounting part [1]	
for boar dummy 11"	11100/0204
for boar dummy 9"	11100/0203
Replacement cover [2]	
for boar dummy 11"	5011100/0233
for boar dummy 9"	5011100/0133
Mobile support for boar dummies [3]	11100/0110
Rubber mat for boar dummy, 150 cm x 100 cm x thickness 23 mm [4]	11100/0040
Semen collection cup, insulated, extra wide opening, volume 1 l [5]	11123/1000
Semen collection bag, 210 x 300 mm, volume 1 l [6]	11123/3010
Filter for boar semen, non-woven tissue, non-absorbing, 240 mm diameter, 200/package [7]	11123/3020



Collection dummy with
stainless steel body