



TECHNICAL DATA SHEET

AUTOMATIC DILUTOR/DISPENSER

Automatic dilutor/dispenser for automated and precise dispensing of two volumina. Use of disposable capillary tubes for sample preparation. With touch screen and memory function.

Product features

- **Accuracy and precision:** The positive displacement system provides better than 99% accuracy, independent of a liquid's viscosity, vapor pressure and temperature.
- **Process security:** The dispenser automates the pipetting process. Methods established by a lab manager are stored in the controller and are consistently run by all technicians.
- **Chemically inert:** The fluid path for the dispenser is made of borosilicate glass, PTFE and CTFE, which are resistant to harsh chemicals used in the laboratory. These high quality materials minimize the potential for leaching of metals or other contaminants that can be detected by downstream processes.

Automatic dilutor/dispenser, preset to 5 ml / 50 µl
dual syringe dilutor with handle [REF. : 12421/0504](#)

