



EasyCoder



The innovation in semen labs: printing of straws without ink

Your benefits

- Unmatched precise and clear printing on 0.5 and 0.25 ml straws
- Up to 3 lines with a total of about 200 characters, even on Ministraws
- Printing of barcodes, logos and special characters
- Immediately abrasion-resistant print
- Control via PC by a clear and user-friendly Minitube software
- Sequential numbers can be printed on straws
- Adjustable character fonts and graduations from micro to macro
- Freely eligible number of lines and of printing positions
- Clean and odorless working without ink and without solvent
- Easy handling, only little maintenance required
- Easy and quick conversion between different straw sizes

Product features

- Thermal transfer printer with automatic straw supply
- Output: up to 3600 straws per hour
- Windows based Minitube software with data base for storing animal information (name, number, ...)
- Data transfer by RS-232 serial port or USB adapter (included)
- Individual composition of printing lines possible, different configurations
- Dimensions (L x W x H): 690 x 370 x 310 mm (PC not included)
- Power supply: 110-230 V AC, 50/60 Hz

EasyCoder printer [REF.: 13038/0000](#)

Color ribbon, 200 m/roll, black [REF.: 13038/0010](#)

