



# Minitube ultrasound scanner

## Portable ultrasound scanner for transabdominal pregnancy detection in sows and small ruminants

### Your benefits

- Compact and light weight
- Portable, battery operated
- Robust, electric probe
- Very easy to clean
- Integrated auto-measure function to easily measure back fat thickness

Stable carrying case, electronic 3.5 MHz convex sensor (70° forward), freeze-function, carrying strap, two rechargeable batteries, battery charger, USB cable to transmit photos

Ultrasound scanner Minitube [REF.: 23500/1000](#)



### Accessories

Spare battery set (2 pc.) [REF.: 5023500/1010](#)

ReproJelly, lubricant, 100 ml [REF.: 11907/0100](#)

### Technical details

|                       |                                       |
|-----------------------|---------------------------------------|
| Monitor size:         | 5.0" TFT-LCD                          |
| Display modes:        | B, B+B                                |
| Image gray scale:     | 256 levels                            |
| Measure functions:    | Distance, circumference, area, volume |
| Detection depth:      | up to 240mm                           |
| Battery capacity:     | 3000mAh, 7.4V                         |
| Power consumption:    | 8W                                    |
| Weight (incl. probe): | 854 g                                 |



Pregnant sow, 22 days



Pregnant sow, 29 days



Pregnant sow, 40 days



Pig fetus, 60 days