



Minitube flush media for embryo transfer



Quality and security guaranteed in every drop

BoviFlush and **EquiFlush** are premium media designed for collecting embryos from cows and mares. Each formulation contains powerful buffering agents, BSA (Bovine Serum Albumin) or PVA (Polyvinyl Alcohol), and antibiotics to ensure optimal embryo recovery and protection.

- + GMP-certified production in a dedicated clean room with stringent production and quality control standards
- + Every batch undergoes veterinarian-supervised quality testing to guarantee safety and performance
- + Ready-to-use with no additional additives required
- + Gas-tight packaging ensures an extended shelf life



Designed for performance, packed for your convenience

The right packaging makes all the difference when it comes to embryo flushing. That's why professionals who value ease of use, precise dosing, hygienic handling, and eco-friendly disposal prefer Minitube **flush media bottles**.

With Minitube flush media, you can count on consistent quality, delivered in packaging designed to make your work easier, cleaner, and more efficient – bottle after bottle.



Best choice – Bottles

Transparent bottles with screw cap and two injection ports [1] allow instant visual inspection and secure, simple handling.

- **One-material construction:** fully recyclable, supporting environmentally responsible disposal
- **Convenient handling:** ergonomic shape with a secure closure for easy transport
- **Hanger included*** [3]: when suspended, the bottle enables gravity-driven emptying of the flushing solution
- **Easy connection:** the spike port of the Y-tubing connects securely to the bottle via one injection port
- **Complete emptying:** the second injection port allows insertion of a cannula at the end of the flushing procedure, enabling air entry and full drainage of the solution

Also available – Infusion bags

For situations where bags are preferred, our **infusion bags** [2] have a three-way port configuration and sealed ports for hygienic fluid removal. However, Minitube bottles remain the best choice for maximum convenience and eco-friendliness.



*Hanger for 1 liter bottle, plastic, multiple use



BoviFlush



BoviFlush with BSA

bag, 1 l	19982/6001
bag, 2 l	19982/6002
bottle with screw cap and two injection ports, 1 l, with hanger	19982/6004

BoviFlush with PVA

bag, 1 l	19982/6011
bag, 2 l	19982/6012
bottle with screw cap and two injection ports, 1 l, with hanger	19982/6010

EquiFlush with BSA

bag, 1 l	19982/6201
bag, 2 l	19982/6202
bottle with screw cap and two injection ports, 1 l, with hanger	19982/6204

EquiFlush with PVA

bag, 1 l	19982/6211
bag, 2 l	19982/6212
bottle with screw cap and two injection ports, 1 l, with hanger	19982/6214

EquiFlush



Minitube's liquid media production

- GMP-compliant pharmaceutical production
- Cleanroom classification ISO 8 / Class D
- Bulk media production ensures large, reliable batches
- Automated bottling under laminar flow
- Contact-free dispensing in a fully sterile environment

Automation and sterile conditions guarantee standardized products of the highest quality

