



Technical Data Sheet

Canine vaginal specula

Product description

Vaginal specula for dogs, model Hannover, stainless steel, consisting of:

- Vaginal specula set of 3 different sizes: diameter 8 mm, 10 mm and 15 mm; length 150 mm
- Light source (3.5 V Halogen)
- Handle with charger for light source
- Scratch-resistant rotatable loupe, 1.6x magnification

Application

Used for vaginoscopy in canine to determine the correct moment of breeding and for AI. Also ideal for collecting vaginal samples.

Properties

- Set of specula for different sized females
- Blunt obturator included
- Bright halogen light, fibre optic
- Handle with charger for light source, can be charged directly in the wall outlet

Order information

Set of vaginal specula	REF.: 17500/0120
Halogen light source	REF.: 17500/0130
Handle with charger	REF.: 17500/0125

