



Suzy, porcine AI training model



A new approach to artificial insemination training

Suzy is a sow model for training AI technicians in a safe and convenient environment before they start working with real animals.

The training model consists of a plastic body and silicone internals. It provides a **most realistic gynecological tract (GT)** that can be removed to demonstrate anatomy and the AI process as an initial step of the training session.

Suzy's GT features a cervix, similar to the natural one, for practicing the lock of the AI catheter in the cervix. A sound will confirm success when the introduced catheter reaches the correct position in the cervix („oink“). A replica of the urethra opening trains the correct introduction of the AI catheter avoiding positioning it in the urethra. Inseminating Suzy with a „real“ semen dose is possible, as the GT is provided with a drain. Both, **traditional** and **post-cervical AI** can be practiced. The GT features a clear tube, where the trainer can overview the advance of the inner cannula of the PCAI catheter.

Suzy provides flexibility around timing and location of AI training courses and lowers the risk of injury to animals by unexperienced trainees. Training units can be repeated as often as necessary. All these benefits help to provide **optimal training** and ultimately an improved performance of applicants of AI in sows.

Product features

- Anatomically correct sow model including full reproductive tract and cervix which provides an authentic feeling during AI training
- Modular approach providing different training levels for quick learning success
- Passage of inner cannula of PCAI catheter is very natural through noticeable cervix cushions
- Correct introduction of AI catheter and inner cannula is confirmed by oinking sound
- „Real“, complete insemination process can be performed
- Compact and light-weight training model for easy transport and cleaning



Increase the efficiency of your AI training courses and contribute to the enhancement of animal welfare!



Order information

Suzy, porcine AI training model - with complete genital tract for traditional and post-cervical insemination

REF.: [22400/0100](#)

Spare parts

(1) Outer part with rectum, vulva, and tail

REF.: [22400/0105](#)

(2) Genital tract for traditional and post-cervical insemination with „Oinkmaster“, vaginal tube and cervix

REF.: [22400/0106](#)



Starter-kit includes

100 AI catheters (Ref. 17106/9176) and
2 tubes of ReproJelly (Ref. 11907/0100)

