



Jute transport box for chilled semen



Sustainable, space-saving, and functional packaging solution

Your benefits

- **Sustainable:** Inlays are made of jute, a naturally renewable raw material, and are fully compostable
- **Space-saving:** The Jute box is much smaller than conventional Styrofoam transport boxes; the cardboard box can be stored folded down with the jute inlays separate from the box
- **Simple:** Easy to pack
- **Functional:** Provides excellent thermal insulation and shock absorbance to protect valuable semen during transport

Product features

- For chilled stallion and dog semen
- With 2 tubes (50 ml) and 2 cool packs
- Outer dimensions: 200 x 185 x 200 mm (W x D x H)
- Holds constant temperature for >30h
- Recommended for transport at ambient temperatures between +5°C and +25°C
- Fully compostable Jute inlays can be disposed of either in the organic waste or the compost

Jute transport box

REF.: [17229/0001](#)



- Styrofoam box -

- Extremely bulky, huge piles of waste
- Has to be disposed of in special recycling facilities

+ Jute box +

- + Strongly reduced waste volume
- + Inlays can be disposed of in the compost or organic waste, the box in the paper waste

+ *Improve your CO2 footprint +*

