



Steripette, capillary pipette



For collecting oocyte cumulus complexes and for transferring oocytes/embryos between dishes

(+ Your benefits

- + Rounded tip for gentle manipulation of oocytes and embryos
- + Plastic capillary cannot break - it is flexible and tough
- + No scratching of the bottom of culture dishes as with glass capillaries
- + Integrated Luer adapter for easy connection to commercially available single use syringes - 1 ml syringes are most convenient
- + Mouse embryo assay (MEA) tested
- + LAL tested for endotoxicity
- + Tested for bioburden

(★ Product features

- Capillary length: 40 mm
- Total length: 60 mm
- Outer diameter of capillary: 1.4 mm
- Inner diameter of capillary: 1.0 mm
- Capillary volume: 40 μ l
- ETO sterilized
- Individually packed for single use

Steripette for manipulation of embryos and oocytes, sterilized, 50/package **19025/0050**

